

PRODUCT DATASHEET

Model: **T-CO**

Description: **T-CO Wall-Mount Transformer Holder**

Revision: January 2020

Dimensions

Product Dimension: 70" x 23.5" x 7.5" (177.80 cm x 59.69 cm x 19.05 cm)

Weight: 66lbs. (29.93 kg)

Boxed Dimension: 75" x 26" x 13" (190.5 x 66.04 x 33.02 cm)

Shipping Weight: 85.5 lbs. (38.78 kg)

FLAMMABILITY RATING

The product is non-flammable under the normal operating conditions

BIOCOMPATIBILITY

This model does not require biocompatible compliance per ISO 10993.

STATEMENT OF CONFORMANCE

Reina Imaging hereby states that the data, specifications, and descriptions contained herein this document are accurate and representative of the models being described.

REGISTRATIONS AND COMPLIANCE

All models listed above are made in the USA and registered to:

Reina Imaging

6107 Lou Ave.

Crystal Lake, IL 60014 USA

AUTHORIZED EU REPRESENTATIVE

PINCHET RX. S.r.l.

20148Milano Via don Gnocchi 24, Italy

+39 02 48706799

DEVICE FAMILY CLASSIFICATION



Country	Class
USA	FDA Class 1, 510(k) exempt
Canada	Class 1
Europe	Class 1, Annex VII, CE
TGA	Class 1